



Dr. perfect V-line



#2-204 Medical Industry Complex, 42-10 Taejunggongdan-gil, Wonju-city, Gangwon-do, 220-962, Republic of Korea
Tel +82-70-8766 2569 Fax +82-33-734-1863 E-mail sales@neodr.com



Dr. perfect V-line





Dr. Perfect V line

Treatment by using absorbable threads that has obtained top grade of CE and KFDA

Possible to perform treatment with application of anesthetic cream

The effect can last for 1~2 years thanks to the induction of collagen synthesis (there may be differences among individuals)

Treatment takes 15minutes~1hour

Able to go back to daily life right after treatment



Dr. perfect V-line


Contents

- 04 INTRODUCTION
- 05 MULTI THREAD
- 06 MONO
- 07 DOUBLE
- 08 TRIPLE
- 09 SCREW(Normal Pitch, Tight Pitch)
- 10 DOUBLE SCREW(Normal Pitch, Tight Pitch)
- 11 COG(COMBI, UNI)
- 12 PACKAGE



What is PDO Therapy

Skin rejuvenation came from the natural wound healing process and PDO therapy originated from traditional thread-lift philosophy in western medicine as well as "threads -embedding therapy" in oriental acupuncture therapy to control and improve tissue stimulation and muscle contraction. The product used is a fine needle (19G-30G) affixed to the absorbable thread (thickness of 0.05mm - 0.40mm and length of 25mm - 140mm). The needle is inserted into the subcutaneous layer of the skin and when it is taken out, the threads remain implanted in the skin. Insertion of dozens of needles to the deeper layer of skin with absorbable threads activates their natural regenerative process and mechanical stimulation of wound healing and tissue contraction



Advantages compared to previous PDO thread Therapy

- 100% PDO (Polydioxanone).
- Clinically proven bio-compatible material that is already widely used in medical suturing.
- Our threads Give both a tissue lift as well as skin rejuvenation.
- Biodegradeable within 6 months and effect can last up to 3 years.

Which areas are treatable using PDO Threads

Upper and Lower Cheeks, Forehead (Brow Lifting), Jawline, Neck, Nasolabial Folds, Marionette Lines, Arms, Thighs and Abdomen

Facial Contouring

PDO Threads can very effectively Strengthen the Contour of the Face, Reduce Roundness, Lift Sagging Skin in the Cheeks, Jowls, Neck and also produce an excellent Eye Brow Lift

Skin Rejuvenation

Enhances Skin Tone and Tightens Pores

Procedure

Local anesthesia is injected into the treatment area. The PDO thread will then be introduced using a cannula or micro-needle into the required area. The procedure takes around 45-60 minutes depending on what treatment the patient requires. The patient is able to fully resume their everyday routine and social commitments within a couple of days. PDO Therapy is a tailor made treatment which can be adjusted to the patients requirements at any time.

Dr. Perfect V line


Is an excellent and effective cosmetic technique to lift and tighten sagging skin tissues. It re-defines facial contours and induces collagen production. The thread forms part of an integrated support structure for the tissue of the face due to collagen synthesis. The lifting effect is at its peak 6 months later when collagen has formed around the threads. The lifting effect lasts between 2-3 years

Only Dr. Perfect V line

Has Multi Thread.
10 strands of thick sutures tied to 4 guiding threads in 23Gage are inserted simultaneously, reducing the treatment time and pain remarkably and thereby maximizing the treatment effect especially for the deep wrinkle such as Nasolabial, Glabella, Forehead wrinkles and even neck folds.

Only Dr. Perfect V line

- Has caps opener and station
- You can open five caps of needle at one try
 - Put them into the station
 - Pick them up one by one to use



MULTI THREAD

Needle and Cannula (Single-use) Type



23G

Needle Gage	Needle Length	Thread Length	USP
23G	38 / 60	25 / 45	5-0 x 14



30G 29G 27G 26G

Needle Gage	Needle Length	Thread Length	USP
30G	25 / 38	30 / 50	7-0
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	6-0
27G	38 / 50 / 60	50 / 70 / 90	5-0
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	5-0




MONO

Needle Type


MONO

Cannula Type (Single-Use) for Eye Treatment



30G

Needle Gage	Needle Length	Thread Length	USP
30G	25	30	7-0



Needle Gage	Needle Length	Thread Length	USP
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	(7-0) x 2
27G	38 / 50 / 60	50 / 70 / 90	(7-0) x 2
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	(7-0) + (6-0)
25G	60 / 90	90 / 150	(6-0) x 2



DOUBLE (D)

Twisted with Two Sutures




TRIPLE (T)

Twisted with Three Sutures




Needle Gage	Needle Length	Thread Length	USP
26G	60 / 90	90 / 150	(7-0) x 3
25G	60 / 90	90 / 150	(7-0) x 2 + (6-0)
23G	60 / 90	90 / 150	(6-0) x 3





Needle Gage	Needle Length	Thread Length	USP
30G	25 / 38	30 / 50	7-0
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	6-0
27G	38 / 50 / 60	50 / 70 / 90	5-0
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	5-0




Normal Type




Tight Pitch Type

SCREW (S)

DOUBLE SCREW (DS)



Needle Gage	Needle Length	Thread Length	USP
29G	25 / 38 / 50 / 60	30 / 50 / 70 / 90	(7-0) X 2
27G	38 / 50 / 60	50 / 70 / 90	(7-0) x 2
26G	38 / 50 / 60 / 90	50 / 70 / 90 / 150	(7-0) + (6-0)
25G	60 / 90	90 / 150	(6-0) x 2



Normal Type



Tight Pitch Type




23G


Needle Gauge	Needle Length	Thread Length	USP
23G	60 / 90	90 / 150	3-0




SHNEO



SNEO




HPNC



PNC

COG (C)

COMBI



COG (C)

UNI



23G



21G




19G

Needle Gauge	Needle Length	Thread Length	USP
23G	60 / 90	105 / 210	(3-0)
21G	60 / 90	140 / 170	(2-0)
19G	100	150	(0-0)




HKSS



HNEO



ORG



SPK





Needle Gauge	Needle Length	Thread Length	USP
21G	60 / 90	140 / 170	2-0
19G	100	150	0-0



COG (C)

UNI CANNULA TYPE




COG (C)

UNI



Needle Gauge	Needle Length	Thread Length	USP
21G	60 / 90	140 / 170	2-0
19G	100	150	0-0






Needle Gage	Needle Length	Thread Length	USP
21G	60 / 90	140 / 170	2-0
19G	100	150	0-0



SPIRAL 3D



SPIRAL 4D



Dr. perfect V-line

PACKAGE



COG & MULTI THREAD
100EA / Individual box



MONO, DOUBLE, TRIPLE,
SCREW, DOUBLE SCREW
200EA / Individual box



CARTON
2,400EA / CTN

COG (C)

UNI CANNULA TYPE

